

Automatic Electric Leg Bag Drainage Device

Freedom Flow 7 Almond Tree Lane Irvine, CA 92612



Concept Development Group

•Freedom Flow is self-contained.

Freedom Flow automatic leg bag drainage system is a completely self-contained unit that easily slips onto the tubing at the bottom of a leg bag. With just the push of a button, *Freedom Flow* opens to allow urine to drain, and then re-clamps the tubing. Urine never touches the *Freedom Flow* device

•Freedom Flow is convenient.

Powered by its own small rechargeable battery, it's perfect for use with all manual wheelchairs or electric wheelchairs. It does not hook up to any outside battery. It slips on and off with ease. It's not a problem during transfers.

•*Freedom Flow* provides independence. Empty your leg bag discreetly and independently. No scheduling headaches and waiting for assistance. At home or on the job, *Freedom Flow* gives you the independence you need.

Freedom Flow promotes your health. Drink all the water/liquids you should be drinking. Drain your bag as often as needed. No more self-imposed fluid limiting. *Freedom Flow* uses a new, ultra-reliable technology to give you convenient, reliable service. *Freedom Flow* was designed by a professor of Industrial Engineering at a major university working in conjunction with a large Rehabilitation center. *Freedom Flow* is designed to be easy to use, and to give you years of reliable, trouble-free service.

The *Freedom Flow* system consists of three main parts:

- 1. A small rectangular unit containing the clamp mechanism and its own internal rechargeable battery;
- 2. A switch;
- 3. A battery charger.

If you require a different type of switch, such as a puff switch, it may be easily replaced or supplied.

Freedom Flow is a superior product with no solenoid or motor to fail. No extra, expensive battery packs are required. There is no cumbersome battery wiring to mess with. You get more value by far from *Freedom Flow*!

The small unit easily slips onto the "oxygen tubing" at the bottom of your leg bag. Just slip it on, locate the operating switch at a convenient location on your wheelchair and lead the single wire cable from the switch to the unit. *Freedom Flow* is ready to use!

Remember, *Freedom Flow* does **not** need any outside power source (no electric wheelchair batteries are needed). Just plug the *Freedom Flow* unit into the included battery charger once every week or two and the rechargeable battery inside *Freedom Flow* will give years of service, more quality and value. *Freedom Flow* comes complete with everything, including its internal battery and battery charger. It is ready to use. Full operating instructions and a warranty are also provided.

For further information or order your *Freedom Flow* today, contact us by mail, phone or email.

To order by phone, call Sue Kamal, Ph.D. at (949) 464-2082.

To order by mail, send this form to: Freedom Flow Concept Development Group 7 Almond Tree Lane Irvine, CA 92612

Or by email at: sue@conceptdg.com

Visit us on the web at: www.conceptdg.com/ff

Please send me	_ Freedom Flow Unit(s)
Total amount enclosed:	<u>_</u> \$
Credit card type: OVisa	OMasterCard (Check)
Credit card #:	
Expiration date:	
Name as it appears on th	e card:

Billing/Shipping address:

City/State/Zip	
Phone	
Email address	